

Enceinte de photostabilité pharmaceutique MSTC-TPS

Réf. : MSTC-150TPS1 • MSTC-500TPS2 • MSTC-500TPS3 • MSTC-1000TPS3 • MSTC-150LTPS • MSTC-500LTPS2



SPÉCIFICATIONS PAR MODÈLE

Modèle	Capacity(L)	Temperature Range	Temperature Uniformity	Temperature Fluctuation	Illuminance	Security Alarm	Shelves(Standard)	Visible Light+UVA	Consumption(W)	Internal Size(W×D×H)
MSTC-150TPS1	150	15~50°C	±2°C(Same level)	±1°C	100~8000Lux, direct setting	Temperature deviation buzzer alarm, independent over-temperature alarm	1		1500	650×490×500
MSTC-500TPS2	500	2~30°C			100~6000Lux, direct setting		2		2000	680×680×1080
MSTC-500TPS3	500						3		2100	1360×490×1360
MSTC-1000TPS3	1000						3		2500	650×490×500
MSTC-150LTPS	150						1		1500	680×680×1080
MSTC-500LTPS2	500						2		2000	

CARACTÉRISTIQUES TECHNIQUES

Enceinte de photostabilité pharmaceutique pour les essais de stabilité à la lumière selon ICH Q1B. Série MSTC-TPS.