


Enceinte de vieillissement UV UATC-UV3

Réf. : UATC-UV3

36 MONTH WARRANTY CE FDA



SPÉCIFICATIONS PAR MODÈLE

Caractéristique	Valeur
Chamber Size (WxDxH)(mm)	1140×650×400
External Size (WxDxH)(mm)	1300×800×1550
Lamp Wavelength (Optional)	UVA-340, Band 315-400 nm
UVB-313,Band 280-315 nm	Irradiance 0.35-1.2W/m
UVB-351, Band 315-400 nm, (dominant peak 351 nm)	Irradiance 0.35-1.2W/m2
Number of Lamps	8
Lamp Power	40W
Lamp Life	More than 1600 hours
Temperature Range	RT + 10°C ~70°C
Humidity Range	75% RH
Temperature Uniformity	± 3°C
Temperature Fluctuation	± 0.5°C
Blackboard Temperature	50°C ~63°C
Distance Between Sample and Lamp	50mm
Center Distance within the Lamp Tube	70mm
Standard Test Piece Size	75×150mm/42 pcs

CARACTÉRISTIQUES TECHNIQUES

Enceinte de vieillissement par rayonnement ultraviolet UATC-UV3, pour simuler le vieillissement des matériaux exposés aux UV.