

Alcoomètre échelle proof selon spécifications Internal Revenue

Réf. : —



SPÉCIFICATIONS PAR MODÈLE

N° de catalogue	Étendue proof	Taille I.R.	Graduation	Tolérance	Longueur approx.
2339	80 to 120%	B	0.5%	±0.5%	230
2340	0 to 20%	F	0.2%	±0.3%	305
2341	20 to 40%	G	0.2%	±0.4%	305
2342	40 to 60%	H	0.2%	±0.4%	305
2343	60 to 80%	I	0.2%	±0.3%	305
2344	75 to 95%	K	0.2%	±0.3%	305
2345	90 to 110%	L	0.2%	±0.3%	305
2346	105 to 125%	M	0.2%	±0.3%	305
2347	125 to 145%	N	0.2%	±0.3%	305
2348	145 to 165%	P	0.2%	±0.3%	305
2349	165 to 185%	Q	0.2%	±0.2%	305
2350	185 to 206%	R	0.2%	±0.2%	305

CARACTÉRISTIQUES TECHNIQUES