

# Alcoomètre échelle proof selon spécifications IRS, avec certificat

Réf. : —



## SPÉCIFICATIONS PAR MODÈLE

| N° de catalogue | Étendue proof | Taille I.R. | Graduation | Tolérance | Longueur approx. |
|-----------------|---------------|-------------|------------|-----------|------------------|
| 2351            | 80 to 120%    | B           | 0.5%       | ±0.5%     | 230              |
| 2352            | 0 to 20%      | F           | 0.2%       | ±0.3%     | 305              |
| 2353            | 20 to 40%     | G           | 0.2%       | ±0.4%     | 305              |
| 2354            | 40 to 60%     | H           | 0.2%       | ±0.4%     | 305              |
| 2355            | 60 to 80%     | I           | 0.2%       | ±0.3%     | 305              |
| 2356            | 75 to 95%     | K           | 0.2%       | ±0.3%     | 305              |
| 2357            | 90 to 110%    | L           | 0.2%       | ±0.3%     | 305              |
| 2358            | 105 to 125%   | M           | 0.2%       | ±0.3%     | 305              |
| 2359            | 125 to 145%   | N           | 0.2%       | ±0.3%     | 305              |
| 2360            | 145 to 165%   | P           | 0.2%       | ±0.3%     | 305              |
| 2361            | 165 to 185%   | Q           | 0.2%       | ±0.2%     | 305              |
| 2362            | 185 to 206%   | R           | 0.2%       | ±0.2%     | 305              |

## CARACTÉRISTIQUES TECHNIQUES