

Densimètre échelle API selon ASTM

Réf. : -



SPÉCIFICATIONS PAR MODÈLE

N° de catalogue	N° ASTM	Étendue API	Graduation	Longueur (mm)
2163A	21H	0 to 6	0.1	165
2164A	22H	5 to 11	0.1	165
2165A	23H	10 to 16	0.1	165
2166A	24H	15 to 21	0.1	165
2167	25H	20 to 26	0.1	165
2168	26H	25 to 31	0.1	165
2169	27H	30 to 36	0.1	165
2170	28H	35 to 41	0.1	165
2171	29H	40 to 46	0.1	165
2172	30H	45 to 51	0.1	165
2173	31H	50 to 56	0.1	165
2174	32H	55 to 61	0.1	165
2175	33H	60 to 66	0.1	165
2176	34H	65 to 71	0.1	165
2177	35H	70 to 76	0.1	165
2178	36H	75 to 81	0.1	165

N° de catalogue	N° ASTM	Étendue API	Graduation	Longueur (mm)
2179	37H	80 to 86	0.1	165
2180	38H	85 to 91	0.1	165
2181	39H	90 to 96	0.1	165

CARACTÉRISTIQUES TECHNIQUES