

Densimètre Brix

Réf. : -



SPÉCIFICATIONS PAR MODÈLE

N° de catalogue	Étendue Brix	Graduation	Longueur approx. (mm)
2196	0 to 12	0.1	330
2197	9 to 21	0.1	330
2198	19 to 31	0.1	330
2199	29 to 41	0.1	330
2200	39 to 51	0.1	330
2201	49 to 61	0.1	330
2202	59 to 71	0.1	330
2203	69 to 81	0.1	330
2204	79 to 91	0.1	330
2205	0 to 35	0.5	330
2206	35 to 70	0.5	330
2207	0 to 70	1.0	330
2208	-5 to 5	0.1	330
2209	5 to 15	0.1	330
2210	15 to 25	0.1	330

