

Réfrigérant à boules (Allihn) avec raccords filetés

Réf. : –



SPÉCIFICATIONS PAR MODÈLE

N° de catalogue	Rodage douille	Rodage cône	Longueur utile (mm)	Raccord	Conditionnement
107/1	–	19/26	200	Screw Thread	1
107/2	–	24/29	200	Screw Thread	1
107/3	–	29/32	200	Screw Thread	1
107/4	–	34/35	200	Screw Thread	1
107/5	–	40/38	200	Screw Thread	1
106/6	–	50/42	200	Screw Thread	1
107/7	–	50/42	250	Screw Thread	1
106/8	19/26	19/26	200	Screw Thread	1
107/9	24/29	24/29	200	Screw Thread	1
107/10	19/26	19/26	300	Screw Thread	1
107/11	24/29	24/29	300	Screw Thread	1
107/12	24/29	24/29	400	Screw Thread	1
107/13	24/29	24/29	500	Screw Thread	1
107/14	34/35	34/35	500	Screw Thread	1
107/15	34/35	34/35	600	Screw Thread	1
107/16	34/35	34/35	750	Screw Thread	1

