

Réfrigérant droit (Liebig)

Réf. : -



SPÉCIFICATIONS PAR MODÈLE

N° de catalogue	Rodage douille	Rodage cône	Longueur utile (mm)	Raccord	Conditionnement
108/1	14/23	14/23	150	Screw Thread	1
108/2	19/26	19/26	200	Screw Thread	1
108/3	24/29	24/29	200	Screw Thread	1
108/4	19/26	19/26	300	Screw Thread	1
108/5	24/29	24/29	300	Screw Thread	1
108/6	19/26	19/26	350	Screw Thread	1
108/7	24/29	24/29	350	Screw Thread	1
108/8	19/26	19/26	400	Screw Thread	1
108/9	24/29	24/29	400	Screw Thread	1
108/10	24/29	24/29	500	Screw Thread	1
108/11	24/29	24/29	600	Screw Thread	1
108/12	24/29	24/29	750	Screw Thread	1
108/13	34/35	34/35	500	Screw Thread	1
108/14	34/35	34/35	600	Screw Thread	1
108/15	34/35	34/35	750	Screw Thread	1
108/16	34/35	34/35	900	Screw Thread	1

N° de catalogue	Rodage douille	Rodage cône	Longueur utile (mm)	Raccord	Conditionnement
108/17	–	19/26	300	Screw Thread	1
108/18	–	24/29	350	Screw Thread	1
108/19	–	34/35	250	Screw Thread	1
108/20	29/32	29/32	250	Screw Thread	1
108/21	29/32	29/32	400	Screw Thread	1

CARACTÉRISTIQUES TECHNIQUES