

Réfrigérant à air

Réf. : –



SPÉCIFICATIONS PAR MODÈLE

N° de catalogue	Longueur utile (mm)	Douille	Cône	Ø extérieur approx. (mm)	Conditionnement
101A/1	200	14/23	14/23	25	1
101A/2	200	19/26	19/26	25	1
101A/3	250	29/29	24/29	25	1
101A/4	400	19/26	19/26	25	1
101A/5	500	24/29	24/29	25	1
101A/6	500	34/35	34/35	25	1
101A/7	600	19/26	19/26	25	1
101A/8	750	24/29	24/29	25	1
101A/9	750	–	19/26	25	1
101A/10	750	–	24/29	25	1
101A/11	1000	–	19/26	25	1
101A/12	1000	–	24/29	25	1

CARACTÉRISTIQUES TECHNIQUES