

Kit densité pétrole L-50, avec certificat NABL

Réf. : -



SPÉCIFICATIONS PAR MODÈLE

N° de catalogue	Composition	Conditionnement
2002L/1	Hydrometer ISI Marked L – 50 SP 0.700 – 0.750 g/ml	1
2002L/2	Hydrometer ISI Marked L – 50 SP 0.750 – 0.800 g/ml	1
2002L/3	Hydrometer ISI Marked L – 50 SP 0.800 – 0.850 g/ml	1
2002L/4	Thermometer Range -10°C to 50°C (Least count. 0.25)	1
2002L/5	OMSONS Filter paper No. 2 (90 mm)	1
2002L/6	Water Paste 100 gms	1
2002L/7	Oil Paste (Dip Paste) 100 gms	1
2002L/8	Glass Jar Borosilicate 500ml.	1

CARACTÉRISTIQUES TECHNIQUES