

Ballon à ébullition à fond plat, col simple à rodage

Réf. : -



SPÉCIFICATIONS PAR MODÈLE

N° de catalogue	Capacité (ml)	Douille	Conditionnement
116/1	5	10/19	10
116/3	10	10/19	10
116/5	25	14/23	10
116/8	50	14/23	10
116/9	50	19/26	10
116/10	50	24/29	10
116/12	100	14/23	10
116/13	100	19/26	10
116/14	100	24/29	10
116/15	100	34/35	2
116/16	150	19/26	10
116/17	150	24/29	10
116/18	150	29/32	2
116/19	150	34/35	2
116/20	250	19/26	10
116/21	250	24/29	10

N° de catalogue	Capacité (ml)	Douille	Conditionnement
116/22	250	29/32	2
116/23	250	34/35	2
116/24	500	19/26	2
116/25	500	24/29	2
116/26	500	29/32	2
116/27	500	34/35	2
116/28	1000	24/29	2
116/29	1000	29/32	2
116/30	1000	34/35	2
116/31	2000	24/29	2
116/32	2000	29/32	2
116/33	2000	34/35	2
116/34	3000	24/29	1
116/35	3000	34/35	1
116/36	3000	40/38	1
116/37	5000	24/29	1
116/38	5000	34/35	1
116/39	5000	40/38	1
116/40	5000	45/40	1
116/41	5000	50/42	1
116/42	10000	34/35	1
116/43	10000	40/38	1
116/44	10000	45/40	1
116/45	10000	50/42	1
116/46	10000	55/44	1
116/47	20000	40/38	1
116/48	20000	50/42	1
116/49	20000	55/44	1

CARACTÉRISTIQUES TECHNIQUES