

Ballon à ébullition à fond rond, col simple à rodage

Réf. : -



SPÉCIFICATIONS PAR MODÈLE

N° de catalogue	Capacité (ml)	Douille	Conditionnement
115/1	5	10/19	10
115/3	10	10/19	10
115/5	25	14/23	10
115/8	50	14/23	10
115/9	50	19/26	10
115/10	50	24/29	10
115/12	100	14/23	10
115/13	100	19/26	10
115/14	100	24/29	10
115/15	100	34/35	2
115/16	150	19/26	10
115/17	150	24/29	10
115/18	150	29/32	2
115/19	150	34/35	2
115/20	250	19/26	10
115/21	250	24/29	10

N° de catalogue	Capacité (ml)	Douille	Conditionnement
115/22	250	29/32	2
115/23	250	34/35	2
115/24	500	19/26	2
115/25	500	24/29	2
115/26	500	29/32	2
115/27	500	34/35	2
115/28	1000	24/29	2
115/29	1000	29/32	2
115/30	1000	34/35	2
115/31	2000	24/29	2
115/32	2000	29/32	2
115/33	2000	34/35	2
115/34	3000	24/29	1
115/35	3000	34/35	1
115/36	3000	40/38	1
115/37	5000	24/29	1
115/38	5000	34/35	1
115/39	5000	40/38	1
115/40	5000	45/40	1
115/41	5000	50/42	1
115/42	10000	34/35	1
115/43	10000	40/38	1
115/44	10000	45/40	1
115/45	10000	50/42	1
115/46	10000	55/44	1
115/47	20000	40/38	1
115/48	20000	50/42	1
115/49	20000	55/44	1

CARACTÉRISTIQUES TECHNIQUES