

Thermomètre usage général

Réf. : -



SPÉCIFICATIONS PAR MODÈLE

N° de catalogue	Étendue	Résolution	Longueur totale	Conditionnement
3004	-10°C to 50°C Mercury Filled	0.50 (1/2)	300	10
3020	-10°C to 50°C Mercury Free Blue Liquid	0.50 (1/2)	300	10
3010	-10°C to 50°C Alcohol Filled	0.50 (1/2)	300	10
3005	-10°C to 110°C Mercury Filled	1.00 (1/1)	300	10
3006	-10°C to 110°C Mercury Filled	0.50 (1/2)	300	10
3021	-10°C to 110°C Mercury Free Blue Liquid	1.00 (1/1)	300	10
3022	-10°C to 110°C Mercury Free Blue Liquid	0.50 (1/2)	300	10
3011	-10°C to 110°C Alcohol Filled	1.00 (1/1)	300	10
3012	-10°C to 110°C Alcohol Filled	0.50 (1/2)	300	10
3007	-10°C to 150°C Mercury Filled	1.00 (1/1)	300	10
3023	-10°C to 150°C Mercury Free Blue Liquid	1.00 (1/1)	300	10
3013	-10°C to 250°C Mercury Filled	1.00 (1/1)	300	10
3008	-10°C to 250°C Mercury Free Blue Liquid	1.00 (1/1)	300	10
3024	-10°C to 360°C Mercury Filled	1.00 (1/1)	300	10
3009	-10°C to 360°C Mercury Filled	1.00 (1/1)	300	10
3009/1	0°F to 120°F Mercury Filled	2.00 (2/1)	300	10

N° de catalogue	Étendue	Résolution	Longueur totale	Conditionnement
3009/2	0°F to 120°F Mercury Free Blue Liquid	1.00 (1/1)	300	10
3019/1	0°F to 120°F Alcohol Filled	1.00 (1/1)	300	10
3009/3	0°F to 220°F Mercury Filled	1.00 (1/1)	300	10
3009/4	0°F to 220°F Mercury Free Blue Liquid	1.00 (1/1)	300	10
3019/2	0°F to 220°F Alcohol Filled	1.00 (1/1)	300	10
3009/5	0°F to 240°F Mercury Filled	1.00 (1/1)	300	10
3009/6	0°F to 240°F Mercury Free Blue Liquid	1.00 (1/1)	300	10
3019/3	0°F to 240°F Alcohol Filled	1.00 (1/1)	300	10
3009/7	0°F to 400°F Mercury Filled	1.00 (1/1)	300	10
3009/8	0°F to 400°F Mercury Free Blue Liquid	1.00 (1/1)	300	10
3019	0°F to 400°F Alcohol Filled	1.00 (1/1)	300	10

CARACTÉRISTIQUES TECHNIQUES