


Éprouvette graduée, embase hexagonale, classe B

Réf. : -



SPÉCIFICATIONS PAR MODÈLE

N° de catalogue	Capacité (ml)	Graduation (ml)	Tolérance (ml)	Bouchon	Conditionnement
80/1	5	0.1	0.1	10/19	2
80/2	10	0.2	0.2	10/19	2
80/3	25	0.5	0.5	14/23	2
80/4	50	1.0	1.0	14/23	2
80/5	100	1.0	1.0	19/26	2
80/6	250	2.0	2.0	24/29	2
80/7	500	5.0	5.0	24/29	1
80/8	1000	10.0	10.0	34/35	1
80/9	2000	20.0	20.0	34/35	1

CARACTÉRISTIQUES TECHNIQUES