

Micropipette multicanaux 8 et 12 canaux

Réf. : –



SPECIFICATIONS PAR MODÈLE

N° de catalogue	Volume (ml)	Incrément ml	Exactitude %	Exactitude ml	CV %	CV ml
MVM-8/10	0.5-10 ml	0.02	1.6	0.16	1	0.1
MVM-8/20	02-20ml	0.02	0.8	0.16	0.4	0.08
MVM-8/50	5-50ml	0.1	0.8	0.4	0.4	0.2
MVM-8/100	10-100ml	0.2	0.8	0.8	0.3	0.3
MVM-8/200	20-200ml	0.2	0.8	1.6	0.3	0.6
MVM-8/300	40-300ml	0.2	0.8	2.4	0.3	0.9
MVM-12/10	0.5-10 ml	0.02	1.6	0.16	1	0.1
MVM-12/20	02-20ml	0.02	0.8	0.16	0.4	0.08
MVM-12/50	5-50ml	0.1	0.8	0.4	0.4	0.2
MVM-12/100	10-100ml	0.2	0.8	0.8	0.3	0.3
MVM-12/200	20-200ml	0.2	0.8	1.6	0.3	0.6
MVM-12/300	40-300ml	0.2	0.8	2.4	0.3	0.9

CARACTÉRISTIQUES TECHNIQUES