

Capsule filtrante PuriCap™ NY-PS (nylon)

Réf. : –

photo

SPÉCIFICATIONS PAR MODÈLE

N° de catalogue	Dimensions (pouce)	Seuil de filtration (µm)	Préfiltre (µm)	Surface filtrante (cm ²)	Conditionnement
OM551-02-02	2	0.20	–	900	1
OM551-05-02	5	0.20	–	1800	1
OM551-08-02	8	0.20	–	2700	1
OM551-02-045	2	0.45	–	900	1
OM551-05-045	5	0.45	–	1800	1
OM551-08-045	8	0.45	–	2700	1
OM551-02-02-08	2	0.20	0.8	700	1
OM551-02-02-045	2	0.20	0.45	700	1
OM551-05-02-08	5	0.20	0.8	1400	1
OM551-05-02-045	5	0.20	0.45	1400	1
OM551-08-02-08	8	0.20	0.8	2100	1
OM551-08-02-045	8	0.20	0.45	2100	1
OM551-02-045-08	2	0.45	0.8	700	1
OM551-05-045-08	5	0.45	0.8	1400	1
OM551-08-045-08	8	0.45	0.8	2100	1

