

Capsule filtrante PuriCap™ PES-PS

Réf. : –



SPÉCIFICATIONS PAR MODÈLE

N° de catalogue	Dimensions (pouce)	Seuil de filtration (µm)	Préfiltre (µm)	Surface filtrante (cm²)	Conditionnement
OM553-02-01	2	0.10	–	500	1
OM553-05-01	5	0.10	–	1000	1
OM553-08-01	8	0.10	–	2000	1
OM553-02-02	2	0.20	–	500	1
OM553-05-02	5	0.20	–	1000	1
OM553-08-02	8	0.20	–	2000	1
OM553-02-045	2	0.45	–	500	1
OM553-05-045	5	0.45	–	1000	1
OM553-08-045	8	0.45	–	2000	1
OM553-02-01-02	2	0.10	0.20	500	1
OM553-02-01-045	2	0.10	0.45	500	1
OM553-05-01-02	5	0.10	0.20	1000	1
OM553-05-01-045	5	0.10	0.45	1000	1
OM553-08-01-02	8	0.10	0.20	2000	1
OM553-08-01-045	8	0.10	0.45	2000	1
OM553-02-02-08	2	0.20	0.80	500	1

N° de catalogue	Dimensions (pouce)	Seuil de filtration (µm)	Préfiltre (µm)	Surface filtrante (cm²)	Conditionnement
OM553-02-02-045	2	0.20	0.45	500	1
OM553-02-02-065	2	0.20	0.65	500	1
OM553-05-02-08	5	0.20	0.80	1000	1
OM553-05-02-045	5	0.20	0.45	1000	1
OM553-05-02-065	5	0.20	0.65	1000	1
OM553-08-02-08	8	0.20	0.80	2000	1
OM553-08-02-045	8	0.20	0.45	2000	1
OM553-08-02-065	8	0.20	0.65	2000	1
OM553-02-045-08	2	0.45	0.80	500	1
OM553-02-045-065	2	0.45	0.65	500	1
OM553-05-045-08	5	0.45	0.80	1000	1
OM553-05-045-065	5	0.45	0.65	1000	1
OM553-08-045-08	8	0.45	0.80	2000	1
OM553-08-045-065	8	0.45	0.60	2000	1

CARACTÉRISTIQUES TECHNIQUES