

Capsule filtrante PuriCap™ PP hydrophile

Réf. : –



SPÉCIFICATIONS PAR MODÈLE

N° de catalogue	Dimensions (pouce)	Seuil de filtration (µm)	Préfiltre (µm)	Surface filtrante (cm ²)	Conditionnement
OM562-02-05	2	0.5	–	500	1
OM562-05-05	5	0.5	–	800	1
OM562-08-05	8	0.5	–	1600	1
OM562-02-10	2	1.0	–	500	1
OM562-05-10	5	1.0	–	800	1
OM562-08-10	8	1.0	–	1600	1
OM562-02-25	2	2.5	–	500	1
OM562-05-25	5	2.5	–	800	1
OM562-08-25	8	2.5	–	1600	1
OM562-02-50	2	5.0	–	500	1
OM562-05-50	5	5.0	–	800	1
OM562-08-50	8	5.0	–	1600	1
OM562-02-100	2	10	–	500	1
OM562-05-100	5	10	–	800	1
OM562-08-100	8	10	–	1600	1
OM562-02-200	2	20	–	500	1

N° de catalogue	Dimensions (pouce)	Seuil de filtration (µm)	Préfiltre (µm)	Surface filtrante (cm²)	Conditionnement
OM562-05-200	5	20	–	800	1
OM562-08-200	8	20	–	1600	1

CARACTÉRISTIQUES TECHNIQUES