

Capsule filtrante PuriCap™ PVDF-PS hydrophile

Réf. : –



SPÉCIFICATIONS PAR MODÈLE

N° de catalogue	Dimensions (pouce)	Seuil de filtration (µm)	Préfiltre (µm)	Surface filtrante (cm²)	Conditionnement
OM557-02-01	2	0.10	–	500	1
OM557-05-01	5	0.10	–	1000	1
OM557-08-01	8	0.10	–	2000	1
OM557-02-02	2	0.20	–	500	1
OM557-05-02	5	0.20	–	1000	1
OM557-08-02	8	0.20	–	2000	1
OM557-02-045	2	0.45	–	500	1
OM557-05-045	5	0.45	–	1000	1
OM557-08-045	8	0.45	–	2000	1
OM557-02-01-02	2	0.10	0.20	500	1
OM557-02-01-045	2	0.10	0.45	500	1
OM557-05-01-02	5	0.10	0.20	1000	1
OM557-05-01-045	5	0.10	0.45	1000	1
OM557-08-01-02	8	0.10	0.20	2000	1
OM557-08-01-045	8	0.10	0.45	2000	1
OM557-02-02-08	2	0.20	0.80	500	1

N° de catalogue	Dimensions (pouce)	Seuil de filtration (µm)	Préfiltre (µm)	Surface filtrante (cm²)	Conditionnement
OM557-02-02-045	2	0.20	0.45	500	1
OM557-05-02-08	5	0.20	0.80	1000	1
OM557-05-02-045	5	0.20	0.45	1000	1
OM557-08-02-08	8	0.20	0.80	2000	1
OM557-08-02-045	8	0.20	0.45	2000	1
OM557-02-045-08	2	0.45	0.80	500	1
OM557-05-045-08	5	0.45	0.80	1000	1
OM557-08-045-08	8	0.45	0.80	2000	1

CARACTÉRISTIQUES TECHNIQUES