

Cartouche filtrante PuriSure™ PES-PS

Réf. : –



SPÉCIFICATIONS PAR MODÈLE

N° de catalogue	Seuil de filtration (µm)	Préfiltre (µm)
OM573-01	0.10	–
OM573-02	0.20	–
OM573-045	0.45	–
OM573-01-02	0.10	0.2
OM573-01-045	0.10	0.45
OM573-01-08	0.10	0.8
OM573-02-045	0.20	0.45
OM573-02-065	0.20	0.65
OM573-02-08	0.20	0.8
OM573-045-65	0.45	0.65
OM573-045-08	0.45	0.8

Type de raccordement	Dimensions (Code) 2.5"(A)	Dimensions (Code) 5"(B)	Dimensions (Code) 10"(C)	Dimensions (Code) 20"(D)	Dimensions (Code) 30"(E)
7P Adaptor (7P)	X	✓	✓	✓	✓
Opti Seal Adaptor (OS)	✓	✓	X	X	X
Seal Clean Adaptor (SC)	✓	✓	X	X	X

N° de catalogue	Adaptateur	Dimensions
OM573-01	7P/OS/SC	A/B/C/D/E

CARACTÉRISTIQUES TECHNIQUES