

Cartouche filtrante PuriSure™ PVDF-PS hydrophile

Réf. : –



SPÉCIFICATIONS PAR MODÈLE

N° de catalogue	Seuil de filtration (µm)	Préfiltre (µm)
OM575-01	0.10	–
OM575-02	0.20	–
OM575-045	0.45	–
OM575-01-02	0.1	0.2
OM575-01-045	0.1	0.45
OM575-02-045	0.2	0.45
OM575-02-08	0.2	0.8
OM575-05-08	0.45	0.8

Type de raccordement	Dimensions (Code) 2.5" (A)	Dimensions (Code) 5" (B)	Dimensions (Code) 10" (C)	Dimensions (Code) 20" (D)	Dimensions (Code) 30" (E)
7P Adaptor (7P)	X	✓	✓	✓	✓
Opti Seal Adaptor (OS)	X	X	X	X	X
Seal Clean Adaptor (SC)	X	X	X	X	X

N° de catalogue	Adaptateur	Dimensions
OM575-02	7P/OS/SC	A/B/C/D/E

CARACTÉRISTIQUES TECHNIQUES