

Ballon à fond rond, quatre cols inclinés

Réf. : -



SPÉCIFICATIONS PAR MODÈLE

N° de catalogue	Capacité (ml)	Douille	Douille latérale	Conditionnement
131/1	100	24/29	14/23	1
131/2	100	24/29	19/26	1
131/3	250	24/29	14/23	1
131/4	250	24/29	19/26	1
131/5	500	24/29	14/23	1
131/6	500	24/29	19/26	1
131/7	500	24/29	24/29	1
131/8	500	34/35	24/29	1
131/9	1000	24/29	14/23	1
131/10	1000	24/29	19/26	1
131/11	1000	34/35	19/26	1
131/12	1000	34/35	24/29	1
131/13	2000	24/29	14/23	1
131/14	2000	24/29	19/26	1
131/15	2000	34/35	19/26	1
131/16	2000	34/35	24/29	1

N° de catalogue	Capacité (ml)	Douille	Douille latérale	Conditionnement
131/17	3000	34/35	19/26	1
131/18	3000	34/35	24/29	1
131/19	5000	34/35	19/26	1
131/20	5000	34/35	24/29	1
131/21	10000	55/44	24/29	1
131/22	20000	55/44	24/29	1

CARACTÉRISTIQUES TECHNIQUES