

Ballon à fond rond, quatre cols parallèles

Réf. : -



SPÉCIFICATIONS PAR MODÈLE

N° de catalogue	Capacité (ml)	Douille	Douille latérale	Conditionnement
132/1	100	24/29	14/23	1
132/2	100	24/29	19/26	1
132/3	250	24/29	14/23	1
132/4	250	24/29	19/26	1
132/5	500	24/29	14/29	1
132/6	500	24/29	19/26	1
132/7	500	24/29	24/29	1
132/8	500	34/35	24/29	1
132/9	1000	24/29	19/26	1
132/10	1000	24/29	19/26	1
132/11	1000	24/29	19/26	1
132/12	1000	34/35	19/26	1
132/13	2000	24/29	14/23	1
132/14	2000	24/29	19/26	1
132/15	2000	34/35	19/26	1
132/16	2000	34/35	24/29	1

N° de catalogue	Capacité (ml)	Douille	Douille latérale	Conditionnement
132/17	3000	34/35	19/26	1
132/18	3000	34/35	24/29	1
132/19	5000	34/35	19/26	1
132/20	5000	34/35	24/29	1
132/21	10000	55/44	24/29	1
132/22	20000	55/44	24/29	1

CARACTÉRISTIQUES TECHNIQUES