

Ballon à fond rond, trois cols parallèles

Réf. : -



SPÉCIFICATIONS PAR MODÈLE

N° de catalogue	Capacité (ml)	Douille	Douille latérale	Conditionnement
126/1	250	24/29	19/26	1
126/2	500	24/29	14/23	1
126/3	500	24/29	19/26	1
126/4	1000	24/29	14/23	1
126/5	1000	24/29	19/26	1
126/6	1000	34/35	24/29	1
126/7	2000	24/29	19/26	1
126/8	2000	34/35	19/26	1
126/9	2000	34/35	24/29	1
126/10	3000	34/35	19/26	1
126/11	3000	34/35	24/29	1
126/12	5000	34/35	19/26	1
126/13	5000	40/38	19/26	1
126/14	5000	50/42	19/26	1
126/15	5000	55/44	19/26	1
126/16	5000	34/35	24/29	1

N° de catalogue	Capacité (ml)	Douille	Douille latérale	Conditionnement
126/17	5000	40/38	24/29	1
126/18	5000	50/42	24/29	1
126/19	5000	55/44	24/29	1
126/20	10000	34/35	19/26	1
126/21	10000	40/38	19/26	1
126/22	10000	50/42	19/26	1
126/23	10000	55/44	19/26	1
126/24	10000	34/35	24/29	1
126/25	10000	50/42	24/29	1
126/26	10000	55/44	24/29	1
126/27	20000	34/35	19/26	1
126/28	20000	50/42	19/26	1
126/29	20000	55/44	19/26	1
126/30	20000	34/35	19/26	1
126/31	20000	34/35	24/29	1
126/32	20000	50/42	24/29	1
126/33	20000	55/44	24/29	1

CARACTÉRISTIQUES TECHNIQUES