

# Ballon à fond rond, deux cols inclinés

Réf. : -



## SPÉCIFICATIONS PAR MODÈLE

| N° de catalogue | Capacité (ml) | Douille | Douille latérale | Conditionnement |
|-----------------|---------------|---------|------------------|-----------------|
| 124/1           | 25            | 19/26   | 14/23            | 1               |
| 124/3           | 50            | 19/26   | 14/23            | 1               |
| 124/5           | 100           | 24/29   | 14/23            | 1               |
| 124/6           | 250           | 24/29   | 14/23            | 1               |
| 124/7           | 250           | 24/29   | 19/26            | 1               |
| 124/8           | 500           | 24/29   | 14/23            | 1               |
| 124/9           | 500           | 24/29   | 19/26            | 1               |
| 124/10          | 1000          | 34/35   | 14/23            | 1               |
| 124/11          | 1000          | 34/35   | 19/26            | 1               |
| 124/12          | 1000          | 24/29   | 14/23            | 1               |
| 124/13          | 1000          | 24/29   | 19/26            | 1               |
| 124/14          | 2000          | 34/35   | 19/26            | 1               |
| 124/15          | 2000          | 34/35   | 24/29            | 1               |
| 124/16          | 2000          | 24/29   | 19/26            | 1               |
| 124/17          | 2000          | 24/29   | 24/29            | 1               |
| 124/18          | 3000          | 34/35   | 24/29            | 1               |

| N° de catalogue | Capacité (ml) | Douille | Douille latérale | Conditionnement |
|-----------------|---------------|---------|------------------|-----------------|
| 124/19          | 3000          | 34/35   | 19/26            | 1               |
| 124/20          | 3000          | 34/35   | 14/23            | 1               |
| 124/21          | 3000          | 24/29   | 24/29            | 1               |
| 124/22          | 3000          | 24/29   | 19/26            | 1               |
| 124/23          | 3000          | 24/29   | 14/23            | 1               |
| 124/24          | 5000          | 40/38   | 24/29            | 1               |
| 124/25          | 5000          | 40/38   | 19/26            | 1               |
| 124/26          | 5000          | 40/38   | 14/23            | 1               |
| 124/27          | 5000          | 34/35   | 24/29            | 1               |
| 124/28          | 5000          | 34/35   | 19/26            | 1               |
| 124/29          | 5000          | 34/35   | 14/23            | 1               |
| 124/30          | 5000          | 24/29   | 24/29            | 1               |
| 124/31          | 5000          | 24/29   | 19/26            | 1               |
| 124/32          | 5000          | 24/29   | 14/23            | 1               |
| 124/33          | 10000         | 55/44   | 34/35            | 1               |
| 124/34          | 10000         | 55/44   | 24/29            | 1               |
| 124/35          | 10000         | 55/44   | 19/26            | 1               |
| 124/36          | 10000         | 50/42   | 34/35            | 1               |
| 124/37          | 10000         | 50/42   | 24/29            | 1               |
| 124/38          | 10000         | 50/42   | 19/26            | 1               |
| 124/39          | 20000         | 55/44   | 34/35            | 1               |
| 124/40          | 20000         | 55/44   | 24/29            | 1               |
| 124/41          | 20000         | 55/44   | 19/26            | 1               |
| 124/42          | 20000         | 50/42   | 34/35            | 1               |
| 124/43          | 20000         | 50/42   | 24/29            | 1               |
| 124/44          | 20000         | 50/42   | 19/26            | 1               |

CARACTÉRISTIQUES TECHNIQUES