

Ballon à fond rond, deux cols parallèles

Réf. : -



SPÉCIFICATIONS PAR MODÈLE

N° de catalogue	Capacité (ml)	Douille	Douille latérale	Conditionnement
125/1	25	19/26	14/23	1
125/3	50	19/26	14/23	1
125/5	100	24/29	14/23	1
125/6	250	24/29	14/23	1
125/7	250	24/29	19/26	1
125/8	500	24/29	14/23	1
125/9	500	24/29	19/26	1
125/10	1000	34/35	14/23	1
125/11	1000	34/35	19/26	1
125/12	1000	24/29	14/23	1
125/13	1000	24/29	19/26	1
125/14	2000	34/35	19/26	1
125/15	2000	34/35	24/29	1
125/16	2000	24/29	19/26	1
125/17	2000	24/29	24/29	1
125/18	3000	34/35	24/29	1

N° de catalogue	Capacité (ml)	Douille	Douille latérale	Conditionnement
125/19	3000	34/35	19/26	1
125/20	3000	34/35	14/23	1
125/21	3000	24/29	24/29	1
125/22	3000	24/29	19/26	1
125/23	3000	24/29	14/23	1
125/24	5000	40/38	24/29	1
125/25	5000	40/38	19/26	1
125/26	5000	40/38	14/23	1
125/27	5000	34/35	24/29	1
125/28	5000	34/35	19/26	1
125/29	5000	34/35	14/23	1
125/30	5000	24/29	24/29	1
125/31	5000	24/29	19/26	1
125/32	5000	24/29	14/23	1
125/33	10000	55/44	34/35	1
125/34	10000	55/44	24/29	1
125/35	10000	55/44	19/26	1
125/36	10000	50/42	34/35	1
125/37	10000	50/42	24/29	1
125/38	10000	50/42	19/26	1
125/39	20000	55/44	34/35	1
125/40	20000	55/44	24/29	1
125/41	20000	55/44	19/26	1
125/42	20000	50/42	34/35	1
125/43	20000	50/42	24/29	1
125/44	20000	50/42	19/26	1

CARACTÉRISTIQUES TECHNIQUES