

# Micropipette monocanal à volume fixe

Réf. : -



## SPÉCIFICATIONS PAR MODÈLE

| N° de catalogue | Volume (ml) | Exactitude $\pm$ % | Exactitude $\pm \mu\text{L}$ | CV $\pm$ % | CV $\pm \mu\text{L}$ |
|-----------------|-------------|--------------------|------------------------------|------------|----------------------|
| FVM-1           | 1.0         | 5                  | 0.05                         | 5          | 0.05                 |
| FVM-2           | 2.0         | 4                  | 0.08                         | 2          | 0.04                 |
| FVM-5           | 5.0         | 2                  | 0.1                          | 1          | 0.05                 |
| FVM-10          | 10.0        | 2                  | 0.1                          | 0.5        | 0.05                 |
| FVM-20          | 20.0        | 0.8                | 0.16                         | 0.4        | 0.08                 |
| FVM-25          | 25.0        | 0.8                | 0.2                          | 0.4        | 0.1                  |
| FVM-50          | 50.0        | 0.8                | 0.4                          | 0.4        | 0.2                  |
| FVM-100         | 100.0       | 0.6                | 0.6                          | 0.2        | 0.2                  |
| FVM-200         | 200.0       | 0.6                | 1.2                          | 0.2        | 0.4                  |
| FVM-250         | 250.0       | 0.6                | 1.5                          | 0.2        | 0.5                  |
| FVM-500         | 500.0       | 0.6                | 3                            | 0.2        | 1                    |
| FVM-1000        | 1000.0      | 0.6                | 6                            | 0.2        | 2                    |
| FVM-2000        | 2000.0      | 0.6                | 12                           | 0.2        | 4                    |
| FVM-5000        | 5000.0      | 0.6                | 30                           | 0.2        | 10                   |
| FVM-10000       | 10000.0     | 0.4                | 40                           | 0.2        | 20                   |

