

Micropipette monocanal à volume variable

Réf. : –



SPÉCIFICATIONS PAR MODÈLE

N° de catalogue	Volume	Incrément $\pm\mu\text{L}$	Exactitude $\pm\%$	Exactitude $\pm\mu\text{L}$	CV $\pm\%$	CV $\pm\mu\text{L}$
VVM-2	0.2-2.0 μL	0.002	2	0.04	1.2	0.024
VVM-10	0.5-10 μL	0.02	1	0.1	0.5	0.05
VVM-20	2-20 μL	0.02	0.8	0.16	0.4	0.08
VVM-50	5-50 μL	0.1	0.8	0.4	0.4	0.2
VVM-100	10-100 μL	0.2	0.6	0.6	0.2	0.2
VVM-200	20-200 μL	0.2	0.6	1.2	0.2	0.4
VVM-1000	100-1000 μL	1.0	0.6	6	0.2	2
VVM-5000	0.5-5 μL	10.0	0.6	30	0.2	10
VVM-10000	1-10 μL	20.0	0.6	60	0.2	20

CARACTÉRISTIQUES TECHNIQUES