

Densimètre à double échelle densité et Baumé

Réf. : -



SPÉCIFICATIONS PAR MODÈLE

N° de catalogue	Étendue de densité	Graduation	Longueur approx.
2279	1.000 to 1.220	0 to 26° Heavy	305
2280	1.000 to 1.225	0 to 26° Heavy	165
2281	1.200 to 1.420	25 to 42° Heavy	305
2282	1.200 to 1.425	42 to 55° Heavy	165
2283	1.400 to 1.620	24 to 43° Heavy	305
2284	1.600 to 1.820	54 to 65° Heavy	305
2285	1.800 to 2.020	64 to 72° Heavy	305
2286	1.000 to 1.600	0 to 45° Heavy	305
2287	1.000 to 1.4	0 to 41° Heavy	165
2288	1.000 to 2.000	0 to 72° Heavy	305
2289	0.700 to 1.000	10 to 100° Light	305

CARACTÉRISTIQUES TECHNIQUES