

Bouchon rodé interchangeable, verre ambré

Réf. : -



SPÉCIFICATIONS PAR MODÈLE

N° de catalogue	Cône	Capacité (mL)	Type
162/1	14/23	Regent Bottle, 60 ml	Flat Head
162/2	19/26	Regent Bottle, 125 & 250 ml	Flat Head
162/3	24/29	Regent Bottle, 500 ml	Flat Head
162/4	29/32	Regent Bottle, 1000 ml	Flat Head
162/5	34/35	Regent Bottle, 2000 ml	Flat Head
162/6	10/19	Volumetric Flask 5,10,20,25 ml	Penny Head
162/7	12/21	Volumetric Flask 50 ml	Penny Head
162/8	14/23	Volumetric Flask 100,200,250 ml	Penny Head
162/9	19/26	Volumetric Flask 500ml	Penny Head
162/10	24/29	Volumetric Flask 1000 ml	Penny Head
162/11	29/32	Volumetric Flask 2000 ml	Penny Head

CARACTÉRISTIQUES TECHNIQUES