


Viscosimètre Cannon-Fenske à écoulement direct

Réf. : –



SPÉCIFICATIONS PAR MODÈLE

N° de catalogue	Dimensions	Nominal Constant	Étendue de viscosité [BS/IP/CF]	Étendue de viscosité [ASTM D 446]	Conditionnement
615/1	25	0.002	0.4 to 1.6	0.5 to 2	1
615/2	50	0.004	0.8 to 3.2	0.8 to 4	1
615/3	75	0.008	NIL	1.6 to 8	1
615/4	100	0.015	3 to 12	3 to 15	1
615/5	150	0.035	7 to 28	7 to 35	1
615/6	200	0.1	20 to 80	20 to 100	1
615/7	300	0.25	50 to 200	50 to 250	1
615/8	350	0.5	100 to 400	100 to 500	1
615/9	400	1.2	240 to 960	240 to 1200	1
615/10	450	2.5	500 to 2000	500 to 2500	1
615/11	500	8	1600 to 6400	1600 to 8000	1
615/12	600	20	4000 to 16000	4000 to 20000	1

CARACTÉRISTIQUES TECHNIQUES