


Viscosimètre en U

Réf. : -



SPÉCIFICATIONS PAR MODÈLE

N° de catalogue	Dimensions	Nominal Constant	Étendue de viscosité [BS/U]	Étendue de viscosité [ASTM D 446]	Conditionnement
612/1	0	0.001	0.3 to 3	NIL	1
612/2	A	0.003	1 to 5	0.9 to 3	1
612/3	B	0.01	2 to 15	2 to 10	1
612/4	C	0.03	6 to 45	6 to 30	1
612/5	D	0.1	20 to 120	20 to 100	1
612/6	E	0.3	60 to 500	60 to 300	1
612/7	F	1	200 to 1000	200 to 1000	1
612/8	G	3	600 to 4000	600 to 3000	1
612/9	H	10	2000 to 15000	2000 to 10000	1

CARACTÉRISTIQUES TECHNIQUES