


Viscosimètre Ubbelohde

Réf. : -



SPÉCIFICATIONS PAR MODÈLE

N° de catalogue	Dimensions	Nominal Constant	Étendue de viscosité [ASTM D 446]	Conditionnement
619/1	0	0.001	0.3 to 1	1
619/2	0C	0.003	0.6 to 3	1
619/3	0B	0.005	1 to 5	1
619/4	1	0.01	2 to 10	1
619/5	1C	0.03	6 to 30	1
619/6	1B	0.05	10 to 50	1
619/7	2	0.1	20 to 100	1
619/8	2C	0.3	60 to 300	1
619/9	2B	0.5	100 to 500	1
619/10	3	1	200 to 1000	1
619/11	3C	3	600 to 3000	1
619/12	3B	5	1000 to 5000	1
619/13	4	10	2000 to 10000	1
619/14	4C	30	6000 to 30000	1

