


Fiole jaugée verre ambré, classe A DIN/ISO

Réf. : -



SPÉCIFICATIONS PAR MODÈLE

N° de catalogue	Capacité (ML)	Col	Tolérance (±ML)	Conditionnement
58/1	1	10/19	0.025	2
58/2	2	10/19	0.025	2
58/3	5	10/19	0.025	1
58/4	10	10/19	0.025	1
58/5	20	10/19	0.04	1
58/6	25	10/19	0.04	1
58/7	50	12/21	0.06	1
58/8	100	14/23	0.1	1
58/9	150	14/23	0.1	1
58/10	200	14/23	0.15	1
58/11	250	14/23	0.15	1
58/12	500	19/26	0.25	1
58/13	1000	24/29	0.4	1
58/14	2000	24/32	0.6	1
58/15	5000	34/35	1.2	1

