

Fiole jaugée verre ambré, classe A, NABL

Réf. : -



SPÉCIFICATIONS PAR MODÈLE

N° de catalogue	Capacité (ML)	Col	Tolérance (±ML)	Conditionnement
59/1	1	10/19	0.025	2
59/2	2	10/19	0.025	2
59/3	5	10/19	0.025	1
59/4	10	10/19	0.025	1
59/5	20	10/19	0.04	1
59/6	25	10/19	0.04	1
59/7	50	12/21	0.06	1
59/8	100	14/23	0.1	1
59/9	150	14/23	0.1	1
59/10	200	14/23	0.15	1
59/11	250	14/23	0.15	1
59/12	500	19/26	0.25	1
59/13	1000	24/29	0.4	1
59/14	2000	24/32	0.6	1
59/15	5000	34/35	1.2	1

