


Fiole jaugée classe A DIN/ISO

Réf. : -



SPÉCIFICATIONS PAR MODÈLE

N° de catalogue	Capacité (ML)	Col	Tolérance ($\pm$ ML)	Conditionnement
50/1	1	10/19	0.025	2
50/2	2	10/19	0.025	2
50/3	5	10/19	0.025	1
50/4	10	10/19	0.025	1
50/5	20	10/19	0.04	1
50/6	25	10/19	0.04	1
50/7	50	12/21	0.06	1
50/8	100	14/23	0.1	1
50/9	150	14/23	0.1	1
50/10	200	14/23	0.15	1
50/11	250	14/23	0.15	1
50/12	500	19/26	0.25	1
50/13	1000	24/29	0.4	1
50/14	2000	29/32	0.6	1
50/15	5000	34/35	1.2	1

