

Fiole jaugée classe A USP, certificat NABL

Réf. : -



SPÉCIFICATIONS PAR MODÈLE

N° de catalogue	Capacité (ML)	Col	Tolérance ($\pm$ ML)	Conditionnement
54/1	5	10/19	0.020	1
54/2	10	10/19	0.020	1
54/3A	20	10/19	0.030	1
54/3	25	10/19	0.030	1
54/4	50	12/21	0.050	1
54/5	100	14/23	0.080	1
54/6	200	14/23	0.100	1
54/7	250	14/23	0.120	1
54/8	500	19/26	0.150	1
54/9	1000	24/29	0.300	1
54/10	2000	29/23	0.500	1

CARACTÉRISTIQUES TECHNIQUES