


Fiole jaugée classe A, NABL

Réf. : -



SPÉCIFICATIONS PAR MODÈLE

N° de catalogue	Capacité (ML)	Col	Tolérance (°ML)	Conditionnement
51/1	1	10/19	0.025	2
51/2	2	10/19	0.025	2
51/3	5	10/19	0.025	2
51/4	10	10/19	0.025	2
51/5	20	10/19	0.04	2
51/6	25	10/19	0.04	2
51/7	50	12/21	0.06	2
51/8	100	14/23	0.1	2
51/9	150	14/23	0.1	2
51/10	200	14/23	0.15	2
51/11	250	14/23	0.15	2
51/12	500	19/26	0.25	2
51/13	1000	24/29	0.4	2
51/14	2000	29/32	0.6	1
51/15	5000	34/35	1.2	1

