

Fiole jaugée classe B

Réf. : -



SPÉCIFICATIONS PAR MODÈLE

N° de catalogue	Capacité (ML)	Col	Tolérance ($\pm$ ML)	Conditionnement
52/1	1	10/19	0.05	2
52/2	2	10/19	0.05	2
52/3	5	10/19	0.05	2
52/4	10	10/19	0.05	2
52/5	20	10/19	0.08	2
52/6	25	10/19	0.08	2
52/7	50	12/21	0.12	2
52/8	100	14/23	0.20	2
52/9	150	14/23	0.20	2
52/10	200	14/23	0.30	2
52/11	250	14/23	0.30	2
52/12	500	19/26	0.50	2
52/13	1000	24/29	0.80	2
52/14	2000	29/32	1.20	1
52/15	5000	34/35	2.50	1

