


Fiole jaugée verre clair, code couleur, classe A

Réf. : -



SPÉCIFICATIONS PAR MODÈLE

N° de catalogue	Capacité (ML)	Col	Tolérance (±ML)	Conditionnement
653/1	1	10/19	0.025	2
653/2	2	10/19	0.025	2
653/3	5	10/19	0.025	1
653/4	10	10/19	0.025	1
653/5	20	10/19	0.04	1
653/6	25	10/19	0.04	1
653/7	50	12/21	0.06	1
653/8	100	14/23	0.1	1
653/9	150	14/23	0.1	1
653/10	200	14/23	0.15	1
653/11	250	14/23	0.15	1
653/12	500	19/26	0.25	1
653/13	1000	24/29	0.4	1
653/14	2000	29/32	0.6	1
653/15	5000	34/35	1.2	1

