

Testeur de point d'éclair Pensky-Martens vase clos FPT-261

Réf. : FPT-261

18 MONTH WARRANTY CE FDA



SPÉCIFICATIONS PAR MODÈLE

Caractéristique	Valeur
Heating Device	The furnace body is made of silicon carbide
Heating Mode	Adjust by manual
Stir Mode	Mechanical drive stirring
Stirring Rate	Procedure A: (90~120) rpm
Oil Cup	Inner diameter: 50.7~50.8 mm
Igniting Source	Gas (or other civilian fuels, the same below)
Thermometers (Mercury-In- Glass Thermometer. Specifications are as Below)	Scale -5~110 °C, division 0.5 °C
Ambient Temperature	15~35 °C
Relative Humidity	≤ 85%
Alimentation	AC (220±10%) V, 50 Hz
Total Power Consumption	≤650 W
Size	370×320×300 mm
N.W.	7 kg
Package Dimension	650×455×360 mm
G.W.	17.6 kg

CARACTÉRISTIQUES TECHNIQUES

Testeur de point d'éclair en vase clos Pensky-Martens FPT-261, pour la détermination du point d'éclair des produits pétroliers, peintures et solvants.