

Appareil de point de fusion MPA-2C / 2A / 1C / 1B

Réf. : MPA-2C • MPA-2A • MPA-1C • MPA-1B



SPÉCIFICATIONS PAR MODÈLE

Modèle	Capillary size	Capillary position capacity	Height of sample	Affichage	Pre-stored measurement mode	Temperature range	Résolution	Linear heating-up rate	Précision	Répétabilité	Linear
MPA-2C	Outer diameter: $\varnothing$ 1.4mm; Inner diameter: $\varnothing$ 1.0mm; Length: 90mm	3	3mm	7-inch color LCD display screen	8	RT ~ 360 °C	0.1 °C	Any continuous options within 0.1°C/min~20°C/min	$\leq$ 200°C $\pm$ 0.4°C; > 200°C $\pm$ 0.7°C	0.3°C (heating rate: 1.0°C/min)	$\pm$ 10%
MPA-2A		1				RT ~ 300 °C					
MPA-1C						RT ~ 400 °C					
MPA-1B						RT ~ 320 °C					

CARACTÉRISTIQUES TECHNIQUES

Appareil de point de fusion MPA-2C / 2A / 1C / 1B, utilisé dans la production de médicaments, colorants, parfums et autres substances cristallines organiques.