

Polarimètre automatique POL-A810 / A850 (P)

Réf. : POL-A810 • POL-A850 • POL-A810P • POL-A850P



SPÉCIFICATIONS PAR MODÈLE

Modèle	Plage de mesure	Measuring Model	Min.Reading	Précision	Repeatability (Standard Deviation)	Light Source	Wavelength	Temperature Control	Temperature Control Range	Précision d
POL - A810	±89.99°Arc; ±259°Z	Optical Rotation; Specific Rotation; Concentration; Sugar Degree	0.001° (optical rotation)	±0.01° (-45°≤optical rotation≤+45°); ±0.02° (optical rotation < -45° or optical rotation > +45°)	0.002°(optical rotation)	LED	589.3nm	\	\	\
POL - A850								Peltier	10°C ~ 50°C	0.3°C
POL - A810P								\	\	\
POL - A850P								Peltier	10°C ~ 50°C	0.3°C

CARACTÉRISTIQUES TECHNIQUES

Polarimètre automatique pour la mesure du pouvoir rotatoire optique (sucres, huiles essentielles, industrie pharmaceutique et alimentaire). Série POL-A810 / A850 / A810P / A850P.