

Testeur de gluten double tête GT-DA

Réf. : GT-DA

18 MONTH WARRANTY CE FDA



SPÉCIFICATIONS PAR MODÈLE

Caractéristique	Valeur
Agitator Rotation Speed	120±2 r/min
Thin Wash Mesh Size	80µm (200mesh) metal sieve
Thick Wash Mesh Size	800µm (26mesh) metal sieve
Washing Liquid Flow Capacity	50mL/min~56mL/min
Dough Mixing Time	1s~99s (digital, adjustable, adjust start from each second)
Washing Time	1s~59s; 1min~15mins (digital, adjustable)
Centrifugation Time	1s~59s; 1min~9mins (digital, adjustable)
Centrifugation Speed	3000 r/min, 6000r/min ± 5r/min (adjustable)
Centrifugation Mesh Box	Aperture is 600µm
Centrifugation Acceleration Speed	2000g
Dryer Temperature	150~ 200 °C
Working Environment	Temperature: 0~40 °C
Alimentation	AC220 V, 50/60 Hz
Shipping Size(W×D×H)(mm)	550×500×460
G.W. (kg)	30

CARACTÉRISTIQUES TECHNIQUES

Testeur de gluten à double tête GT-DA, pour la mesure du gluten humide des farines.