

Testeur de gluten simple tête GT-SA

Réf. : GT-SA



SPÉCIFICATIONS PAR MODÈLE

Caractéristique	Valeur
Agitator Rotation Speed	120±2r/min
Thin Wash Mesh Size	80µm (200 mesh) metal sieve
Thick Wash Mesh Size	800µm (26 mesh) metal sieve
Washing Liquid Flow Capacity	50mL/min~56mL/min
Dough Mixing Time	1s~59s (digital, adjustable, adjust start from each second)
Washing Time	1s~59s; 1min~9mins (digital, adjustable)
Centrifugation Time	1s~59s; 1min~9mins (digital, adjustable)
Centrifugation Speed	3000 r/min, 6000 r/min±5r/min (adjustable)
Centrifugation Mesh Box	Aperture is 500µm
Centrifugation Acceleration Speed	2000g
Dryer Temperature	150°C~ 200°C
Working Environment	Temperature 0°C~40°C
Alimentation	AC220V, 50/60Hz
Shipping Size(W×D×H)(mm)	480×470×440
G.W. (kg)	25

CARACTÉRISTIQUES TECHNIQUES

Testeur de gluten à simple tête GT-SA, pour la mesure du gluten humide des farines.