

Viscosimètre Burgers VSC-B5

Réf. : VSC-B5

12 MONTH WARRANTY CE FDA



SPÉCIFICATIONS PAR MODÈLE

Caractéristique	Valeur
Viscosity measurement range	0-99.99 mPa.s
Viscosity value display resolution	0.01 mPa.s
Instrument measurement accuracy	±0.02 mPa.s
Time measurement range	0-5000 s
Time measurement precision	≤ 0.01 s
Time display resolution	Milliseconds
Temperature control range of water bath	R.T. -70 °C
Water bath temperature control setting resolution	0.1 °C
Temperature control accuracy of water bath	±0.1 °C
Water bath temperature display resolution	0.01 °C
Viscometer volume	100 mL
Volumetric error of the viscometer	± 0.5 mL
Longitudinal center moment of measuring device	152 mm ± 0.5 mm
Alimentation	110/220 V, 50/60 Hz
Consumption	500 W
Shipping Size	465×370×610 mm

CARACTÉRISTIQUES TECHNIQUES

Viscosimètre Burgers VSC-B5 pour la mesure de la viscosité selon la méthode à coupe, adapté aux peintures et revêtements.