

Viscosimètre numérique VSC-1N / 2N / 3N / AIN

Réf. : VSC-1N • VSC-2N • VSC-3N • VSC-AIN

24 MONTH WARRANTY CE FDA



SPÉCIFICATIONS PAR MODÈLE

Modèle	Measuring Range (mPa·s)	Rotating Speed (RPM)	Number of Spindles	Précision	Répétabilité	Readability	Alimentation	LCD Screen Display	L0 Spindle is for Option Measured Range (mPa·s)
VSC-1N	20-600,000	1-60	4 types(L1~L4)	±1%(FS) (Newtonian Fluid)	±0.5%(FS) (Newtonian Fluid)	1. when viscosity is lower than 10.0 cp: 0.01 2. when viscosity is lower than 10000 cp: 0.1 3. when viscosity is equal or higher than 10000 cp: 1	AC100~240V, 50/60Hz	The LCD screen can display the viscosity, speed, torque, spindle and the maximum viscosity can be measured in the current spindle speed value.	1~600
VSC-2N	12-6,000,000	0.1-99.9	4 types (21, 27, 28, 29)						0.6~6,000
VSC-3N	160-80,000,000	0.1-99.9							6~60,000
VSC-AIN	50-10,000,000	0.1-200							3~60,000

CARACTÉRISTIQUES TECHNIQUES

Viscosimètre numérique pour la mesure de la viscosité des liquides en laboratoire et en production. Affichage numérique. Série VSC-1N / 2N / 3N / AIN.