

Viscosimètre rotatif numérique VSC-LDV / RDV / HDV

Réf. : [VSC-LDV1](#) • [VSC-RDV1](#) • [VSC-HDV1](#) • [VSC-LDV2](#) • [VSC-RDV2](#) • [VSC-HDV2](#)



SPÉCIFICATIONS PAR MODÈLE

Modèle	Series	Measuring Range(mPa·s)	Rotating Speed (RPM)	Number of Spindles	Précision	Répétabilité	Readability	Alimentation	LCD Screen Display	L0 Spindle is for Option Measured Rang
VSC-LDV1	DV-1	12~6M	0.1 ~ 99.9	4 types (L1~L4)	±1%(FS) (Newtonian Fluid)	±0.5%(FS) (Newtonian Fluid)	1. when viscosity is lower than 10.0 cp : 0.01 2. when viscosity is lower than 10000 cp : 0.1 3. when viscosity is equal or higher than 10000 cp : 1	AC100~240V, 50/60Hz	The LCD screen can display the viscosity, speed, torque, spindle and the maximum viscosity can be measured in the current spindle speed value.	0.6~6K
VSC-RDV1	DV-2	80~40M	0.1 ~ 200	6 types(R2~R7)						6~60K
VSC-HDV1		640~320M		4 types (L1~L4)						50~500K
VSC-LDV2		6~6M		6 types(R2~R7)						0.3~6K
VSC-RDV2		40~40M								3~60K
VSC-HDV2		320~320M								25~500K

CARACTÉRISTIQUES TECHNIQUES

Viscosimètre rotatif numérique à plages de mesure étendues (basse, moyenne et haute viscosité). Série VSC-LDV / RDV / HDV (versions 1 et 2).